

Determination of artemisinin content by hand-held NIR and FT-NIR spectroscopy: A promising tool for breeding.

C. Camps, P. Rauscher, J. Henkel, X. Simonnet, M. Toussiot, M. Quennoz

¹Agroscope, Centre de recherche conthey, CH-1964 Conthey, Switzerland; ²Hochschule Hannover, University of Applied Sciences and Arts, Germany. cedric.camps@acw.admin.ch

Introduction

Artemisinin is a drug used to treat malaria. It is a sesquiterpene lactone with an endoperoxide bridge which is difficult to synthesize in a cost-effective manner, so direct extraction from the leaves of *A. annua* remains an essential step in drug manufacture. The aim of the present study, models aiming at determining the artemisinin content (AC) were evaluated using a 2 NIR spectrometers: a hand-held (phazir 10-18, polychromix)¹ and a lab unit (FT-NIR, MPA Bruker Optics).

Calibration

A set of 60 samples of *A. annua* dried leaves were used for calibrations. PLS regressions including cross-validation were used to correlate NIR data and references AC-values. Reference AC-values have been determined by thin layer chromatography. The table 1 presents the calibration (with cross-validation) results and the figure 1 shows how the reference and the predicted values fit in the calibration step using the handheld device.

TABLE 1. PLS-values of cross-validation obtained with the two NIR devices.

PLS paramters	Calibration with Cross-validation				
	R ²	SECV	LV	λ range	RPD
Hand-held NIR	0.94	0.1	6	948-1781 nm	2.32
FT-NIR	0.92	0.03	10	7506-6094; 5454.1-4597.8 cm ⁻¹	3.49

R²: coefficient of determination, SECV: standard error of cross-validation, LV: number of latent variables, λ range: the range of wavelenghts used in the model, RPD: ration performance to deviation.




Fig.1. Actual vs. Predicted values of AC (%) using the handheld NIR device (cross-validation).

Validation

Validations were performed for both devices. SEP-values have been calculated without and with bias correction. PLS-values of validation steps are presented in the table 2. Figure 2 shows a measurement performed with the hand-held device and the display of the content of artemisinin.

TABLE 2. PLS-values of validation obtained with the two NIR devices.

PLS paramters	Validation		
	SEP	bias	SEPC
Hand-held NIR	0.103	-0.006	0.103
FT-NIR	0.062	0.02	0.060

SEP: standard error of prediction, bias: bias value, SEPC: SEP corrected from bias.



Fig.2. Actual vs. Predicted values of AC (%) using the FT-NIR device (validation).

Conclusion

- Prediction with the FT-NIR is satisfactory.
- The results obtained with the hand-held NIR spectrometer are promising but not ready for use in practice because the values of RPD and SEP.
- Develop an accurate and robust calibration with the hand-held device requires a larger number of samples and a much larger time labor and work compare to F-NIR device.
- Within the framework of a breeding program aiming at selecting high content hybrids, FT-NIR and hand-held NIR spectroscopy remain interesting methods compared to chemical analyses (TLC, HPLC).



Schweizerische Eidgenossenschaft
Confédération suisse
Confederazione Svizzera
Confederaziun svizra

Federal Department of Economic Affairs,
Education and Research EAER
Agroscope

Swiss Confederation